

Three-field Ionization Chamber Calibration Procedure for
 Pre-Amplifier board assemblies 61170A
 (Also valid for the pre-amp assemblies listed below)

The following adjustments apply to the calibration of a 61170A pre-amplifier board for a stationary 3-field ion chamber, e.g. for chest or table use.

Pre-amp Assembly	Description	Difference from 61170A
61170A	Calculated Gain Range = 1 to 21	None
61170C	Calculated Gain Range = 1.47 to 34.8	R21 = 1500 ohm R25 = 2200 ohm R4 = 50K ohm
61170J	Calculated Gain Range = 1.47 to 34.8	R21 = 1500 ohm R25 = 2200 ohm R4 = 50K ohm Side-Turn Pots
61170N	Calculated Gain Range = 2.2 to 52.2	R21 = 1000 ohm R25 = 2200 ohm R4 = 50K ohm
61170R	Calculated Gain Range = 2.1 to 24.9	R21 = 2200 ohm R25 = 4700 ohm R4 = 50K ohm Side-Turn Pots
61170S	Calculated Gain Range = 2.68 to 63.7	R21 = 820 ohm R25 = 2200 ohm R26 = 680 ohm R4 = 50K ohm
61170T	Calculated Gain Range = 2.68 to 63.7	R21 = 820 ohm R25 = 2200 ohm R26 = 680 ohm R4 = 50K ohm Side-Turn Pots

WARNING:

- Do not attempt to service the equipment unless this service manual has been consulted and is understood.
- Failure to heed this warning may result in injury to the service provider, operator or patient from electric shock and mechanical or other hazards.

Note: When working with the pre-amplifier assembly it is important that electrostatic discharge (ESD) prevention techniques be observed. Before touching the pre-amp assembly, attach an ESD wrist strap to yourself. Be sure to ground yourself and the ion chamber frame to dissipate static charges.

Note: The pre-amp assembly is a very delicate and sensitive device. It is important to keep it as clean as possible. Wash and dry your hands thoroughly before working with it and, when possible, use unpowdered latex or cotton gloves. Take care to touch the pre-amp board as little as possible. Take extra care to avoid touching the three air-mounted field inputs. Oils from your fingers on the air-mounts or their components can cause performance degradation.

This document represents proprietary information originated by Advanced Instrument Development, Inc. and which shall not be disclosed or utilized in any manner detrimental to the company's business.

The procedure assumes that the installation of the Automatic Exposure Control (AEC) is complete and that the AEC and x-ray generator are in proper working condition. After making the necessary interconnections between the ion chamber and the AEC, power up the system.


The 61170A pre-amplifier board can be configured to operate with several different AEC systems. If specified at the time of purchase, the ion chamber will be delivered with the pre-amplifier board already configured for a particular application. To reconfigure an ion chamber for a different configuration, see the section on Ion Chamber Inputs and Output and the section on Specific Configurations at the end of this document.

Calibration/Test Set Up:

Your X-Ray Generator/AEC Manual should cover this subject. Calibration procedures are specific to the X-Ray Generator/AEC. The Following is provided as additional information only.

Select the center field of the ion chamber. Set the generator for 100 kVp and maximum backup time. For 100 kVp use 8 to 10 inches (20 to 25 cm) of water or plastic for a phantom. Metals such as copper, aluminum or lead are not suitable for use as phantoms. Make sure the phantom is homogeneous and completely covers all fields equally. Center the x-ray beam on the center field. Collimate the x-ray beam so that it completely covers all three fields but does not extend beyond the limits of the phantom.

Adjustment Potentiometers:



All necessary adjustment potentiometers are accessible through the pre-amp chassis cover. There is no need to open the pre-amp cover during normal calibration procedures.

Master Gain Adjustment:

The master gain adjustment should be used to adjust the ion chamber to produce the correct optical density (OD) while only the Center Field (Field B) is selected. A replacement ion chamber's film optical density should usually be matched to that of other stationary ion chambers connected to the system. Note that the master gain adjustment is a multi-turn potentiometer. A clockwise adjustment to the master gain potentiometer will increase the sensitivity of the chamber, causing the length of the exposure (mAs) to decrease.

Make an exposure and record the post exposure mAs reading. Process the film. Measure and record the film optical density. Adjust the master gain for the desired optical density (normally between 1.1 and 1.8 OD, depending upon the user). Make the master gain adjustment using this formula:

$$\text{Desired mAs} = (\text{Desired OD} / \text{Measured OD}) \times \text{Measured mAs}$$

Note: Make sure that a film cassette is present while making all exposures. Cassette backscatter affects AEC operation.

Field Balance Check:

Using the AEC post-exposure mAs display or other calibrated mAs meter check the Left Field (Field A) and Right Field (Field C) individually to see that they each produce the same mAs reading as the Center Field (Field B).

If necessary, adjust the individual gain potentiometers to balance the outputs to give the same mAs reading for each field selected individually. Do not adjust Field B unless you run out of gain range on the Left and Right Field adjustments. If Field B requires readjustment, be certain to repeat the Master Gain adjustment.

Note that individual gain adjustments are multi-turn potentiometers. A clockwise adjustment to a field gain potentiometer will increase the sensitivity of that field, causing the length of the exposure (mAs) to decrease.

Ionization Chamber Pin-outs:

61170A Pre-amp Board Pin-out	Function	9-Pin Sub-D Pin-Out
1	NONE	1
2	FIELD 2 SELECT	2
3	FIELD 1 SELECT	3
4	RESET	4
5	OUTPUT	5
6	FIELD 3 SELECT	6
7	NEGATIVE SUPPLY	7
8	POSITIVE SUPPLY	8
9	GROUND	9


NOTE: Cable wire colors do not match those inside the pre-amp chassis. Cable pin-out details are available on-line at <http://www.aidxray.com> or by contacting Advanced Instrument Development, Inc.

Acceptable Power Supply Ranges for 61170A Pre-amp:

Supply Voltage	Measurement Point	Acceptable Range
External +12VDC	61170A pin 8 (Referenced to 61170A pin 9)	+11.4VDC to +15.8VDC
External -12VDC	61170A pin 7 (Referenced to 61170A pin 9)	-11.4VDC to -15.8VDC
Internal +12VDC Isolated	Measure across 61170A C23	+10.8VDC to +12.5VDC
Internal -12VDC Isolated	Measure across 61170A C22	-10.8VDC to -12.5VDC
Internal +5VDC	Measure across 61170A C18	+4.7VDC to +5.3VDC
Internal +75VDC	Measure across 61170A C49 (HV+ to HV-)	+65VDC to +85.0VDC

Ionization Chamber Inputs and Output:

Signal	Switch	Comments
Positive Supply +12VDC	None	+11.4VDC to +15.75VDC less than 85 mA.
Negative Supply -12VDC	None	-11.4VDC to -15.75VDC less than 15 mA.
Low-Active Commands Reset & Field Selects	SW1-1 OFF SW1-2 OFF SW1-3 OFF SW1-4 ON	Pulling the Reset line to ground ($\leq 6VDC$) beginning at exposure start and lasting for the entire duration of the exposure allows the integrator to operate. Pulling the field select lines to ground ($\leq 6VDC$) will select the field.
High-Active Commands Reset & Field Selects	SW1-1 OFF SW1-2 OFF SW1-3 ON SW1-4 OFF	Driving the Reset line high ($\geq 8VDC$) beginning at exposure start and lasting for the entire duration of the exposure allows the integrator to operate. Driving the field select lines high ($\geq 8VDC$) will select the field.
Field Configuration A=1, B=2, C=3	SW2-1,5 & 8 ON SW2-2,3,4, 6 & 7 OFF	Field selects 1, 2 and 3 select left (A), center (B) and right (C), respectively, as viewed from the x-ray tube-side of the ion chamber.
Field Configuration A=2, B=1, C=3	SW2-2,4 & 8 ON SW2-1,3,5, 6 & 7 OFF	Field selects 2, 1 and 3 select left (A), center (B) and right (C), respectively, as viewed from the x-ray tube-side of the ion chamber.
Field Configuration A=3, B=1, C=2	SW2-3,4 & 7 ON SW2-1,2,5, 6 & 8 OFF	Field selects 3, 1 and 2 select left (A), center (B) and right (C), respectively, as viewed from the x-ray tube-side of the ion chamber.
Field Configuration A=3, B=2, C=1	SW2-3,5 & 6 ON SW2-1,2,4, 7 & 8 OFF	Field selects 3, 2 and 1 select left (A), center (B) and right (C), respectively, as viewed from the x-ray tube-side of the ion chamber.
Positive Output	Not Selectable Specified at time of purchase	Time integrated signal ramping up from 0VDC to a maximum of at least 80% of the positive supply voltage (+9.6VDC for +12VDC supply). The slope of this signal is directly proportional to the amount of x-ray flux received.
Negative Output	Not Selectable Specified at time of purchase	Time integrated signal ramping down from 0VDC to a maximum of at least 80% of the negative supply voltage (-9.6VDC for -12VDC supply). The slope of this signal is directly proportional to the amount of x-ray flux received.



Specific Configurations:

Unless specified otherwise, ICX ion chambers are delivered with an AID compatible jumper configuration. AID compatible means that the input and output signals will interface with Advanced Instrument Development, Inc's Expos-AID™ Automatic Exposure Control. This same configuration will also interface with Acoma, Control-X, CPI, Del Medical (Gendex), Electromed (EMD, Triton), OEC, Quantum Medical Imaging, Sedecal (Innerscan), Summit Industries, etc.

AID Compatible jumper configuration:

Function	Switch Setting
Low-Active Commands Reset & Field Selects	SW1-1 OFF SW1-2 OFF SW1-3 OFF SW1-4 ON
Field Configuration A=2, B=1, C=3	SW2-2,4 & 8 ON SW2-1,3,5, 6 & 7 OFF
Positive Output	Not Selectable Specified at time of purchase

Alternate configurations are listed below.

GE, Fischer, Varian Compatible jumper configuration:

Function	Switch Setting
High-Active Commands Reset & Field Selects	SW1-1 OFF SW1-2 OFF SW1-3 ON SW1-4 OFF
Field Configuration A=1, B=2, C=3	SW2-1,5 & 8 ON SW2-2,3,4, 6 & 7 OFF
Positive Output	Not Selectable Specified at time of purchase

Trex Medical - Continental, Keithley, Ratheon, Xonics Compatible jumper configuration:


Function	Switch Setting
Low-Active Commands Reset & Field Selects	SW1-1 OFF SW1-2 OFF SW1-3 OFF SW1-4 ON
Field Configuration A=1, B=2, C=3	SW2-1,5 & 8 ON SW2-2,3,4, 6 & 7 OFF
Positive Output	Not Selectable Specified at time of purchase

GTR Labs Compatible jumper configuration:

Function	Switch Setting
Low-Active Commands Reset & Field Selects	SW1-1 OFF SW1-2 OFF SW1-3 OFF SW1-4 ON
Field Configuration A=3, B=1, C=2 (= GTR Labs fields 1-2-3)	SW2-3,4 & 7 ON SW2-1,2,5, 6 & 8 OFF
Positive Output	Not Selectable Specified at time of purchase

Picker Compatible jumper configuration:

Function	Switch Setting
Low-Active Commands Reset & Field Selects	SW1-1 OFF SW1-2 OFF SW1-3 OFF SW1-4 ON
Field Configuration A=2, B=1, C=3	SW2-2,4 & 8 ON SW2-1,3,5, 6 & 7 OFF
Negative Output	Not Selectable Specified at time of purchase



- NOTES:
1. * MANUFACTURER MAY USE ALTERNATE COMPONENTS UNLESS NOTED OTHERWISE.
 2. ALL SWITCHES OFF UNLESS OTHERWISE NOTED.
 3. U7, U8, U9, AND U10 MAY BE (4) SINGLE, (2) DUAL, OR (1) QUAD DEVICES.
 4. ALL DIODES ARE FMMD914 UNLESS NOTED OTHERWISE

DO NOT SCALE				
UNLESS SPECIFIED OTHERWISE; DIMENSIONS ARE IN INCHES. FRACTIONS DECIMALS ANGLES ±1/64 ±0.010 ±1°				
02	2192	02SEP04	ADDED G3, CONNECTED A5 TO GND	TRW
01	2192	02JUL04	BIAS SUPPLY ISOLATED MULTIPLIER	TRW
00	2192	10MAY04	PROTOTYPE	TRW
REV.	ECN NO.	DATE	REVISIONS	BY

ADVANCED INSTRUMENT DEVELOPMENT, INC.
1011 N. 25TH AVE.
MELROSE PARK, IL. 60160

ISOLATED 3-FIELD PRE-AMP
GAIN RANGE: 1 TO 21
TOP-TURN POTS

SCALE	DRAWN BY	DATE
MATERIAL	CHECKED BY	
FINISH	APPROVED BY	
THIS DRAWING REPRESENTS PROPRIETARY AND CONFIDENTIAL INFORMATION ORIGINATED BY ADVANCED INSTRUMENT DEVELOPMENT, INC. AND WHICH SHALL NOT BE DISCLOSED OR UTILIZED IN ANY MANNER DETRIMENTAL TO THE COMPANY'S BUSINESS.		
USED ON	NEXT ASSY.	DRAWING NO.
ICX SERIES		61170A
REV.		02