

## Section D



# Picker Compatible Ion Chambers

17 January 2011

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

ADVANCED INSTRUMENT DEVELOPMENT, INC. 2545 CURTISS STREET DOWNERS GROVE, ILLINOIS 60515-4059  
TELEPHONE: (630) 271-9729 FACSIMILE: (630) 271-9995 [www.aidxray.com](http://www.aidxray.com)

# Section D



<b>MODEL:</b> ICX143	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154N	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.28 IN (7.1mm)	
<b>APPLICATION:</b> GENDEX / PICKER TABLE BUCKY L.H. LOAD	
AID COMPATIBLE	

<b>MODEL:</b> ICX144	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154N	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.28 IN (7.1mm)	
<b>APPLICATION:</b> GENDEX / PICKER TABLE BUCKY R.H. LOAD	
AID COMPATIBLE	

<b>MODEL:</b> ICX167	
<b>FIELDS:</b> 3-FLD SPLIT CENTER	
<b>PRE-AMP:</b> 61135F	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> 3-FIELD WITH SPLIT CENTER FOR FORMATTED SPOT FILM	
AID COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX183	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154G	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> 14 X 36 CASSETTE R.H. LOAD  PICKER COMPATIBLE	

<b>MODEL:</b> ICX195	
<b>FIELDS:</b> A=2, B=1, (B+C)=3	
<b>PRE-AMP:</b> 61112A	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> PICKER SPOT FILMER  PICKER COMPATIBLE	

<b>MODEL:</b> ICX197	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154N	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.28 IN (7.1mm)	
<b>APPLICATION:</b> PICKER TABLE BUCKY  PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX801	
<b>FIELDS:</b> SPOT (2 IN SQ.)	
<b>PRE-AMP:</b> 61125B	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> EXIT-TYPE PEDIATRIC / SPOT FILM DEVICE	
<b>PICKER COMPATIBLE</b>	

<b>MODEL:</b> ICX811L	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154L	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> POERSCH METALS LH LOAD CHEST STAND	
<b>PICKER COMPATIBLE</b>	

<b>MODEL:</b> ICX818	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61135A	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.28 IN (7.1mm)	
<b>APPLICATION:</b> FUJI FCR9501, FCR9501HQ	
<b>PICKER COMPATIBLE</b>	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX821	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.34 IN (8.4mm)	
<b>APPLICATION:</b> FUJI 9502/HQ TABLE	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX822	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154S	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.28 IN (7.1mm)	
<b>APPLICATION:</b> CANON DIGITAL RADIOGRAPIC SYSTEM	
REMOVABLE REMOTE PRE-AMP PICKER COMPATIBLE	

<b>MODEL:</b> ICX823	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.22 IN (5.7mm)	
<b>APPLICATION:</b> FUJI FCR 9501-ES CHEST STAND	
PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX824	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.22 IN (5.7mm)	
<b>APPLICATION:</b> FUJI FCR 5501 17 X 17 CHEST STAND	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX825	
<b>FIELDS:</b> ICX825L: A=2, B=1, C=3 ICX825R: A=3, B=1, C=2	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.22 IN (5.7mm)	
<b>APPLICATION:</b> FUJI FCR 5502 TABLE	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX826	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61170S	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.13 IN (3.3mm)	
<b>APPLICATION:</b> FUJI VELOCITY 17 X17 CHEST STAND	
PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX846	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154N	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.28 IN (7.1mm)	
<b>APPLICATION:</b> ANALOGIC, KODAK DIGITAL CHEST STAND	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX848	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61116A	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> BIODEX TRAUMA TABLE	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX853R	
<b>FIELDS:</b> A=3, B=1, C=2	
<b>PRE-AMP:</b> 61154J	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> GENERIC R.H. LOAD SIDE -TURN POTS	
PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX854	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154S	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.13 IN (3.3mm)	
<b>APPLICATION:</b> KONICA REGIUS 350 CR SYSTEM	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX855	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> PAUSCH KS80 SIDE-TURN POTS	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX859	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> L-F TABLE WITH GRID MOTOR AT FOOT-END	
PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX862	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> GENERIC L.H. LOAD L-F 9000 R.H. LOAD	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX874	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> L-F HUT III OEC UROVIEW 2600	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX877	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> REPLACEMENT FOR FISCHER PICKER ELITE 1500 (Cat. No. 1280) PICKER ELITE 5000 (Cat. No. 1281)	
PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX890	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61150C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> PAUSCH KS-80	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX892	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> L-F 9000 R.H. LOAD	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX919	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> PICKER AEC, H.L. LYONS R.H. LOAD CHEST STAND	
PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX921	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61135B	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> EXPLORATOR 35 SPOT FILM DEVICE	
AID COMPATIBLE REMOTE	

<b>MODEL:</b> ICX923	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154C	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> PAUSCH KS-80	
PICKER COMPATIBLE	

<b>MODEL:</b> ICX925	
<b>FIELDS:</b> A=1, B=2, C=3	
<b>PRE-AMP:</b> 60917F	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> EXIT-TYPE REPLACES PICKER PMT IN A CHEST STAND FILM CHANGER USING KEITHLEY AEC	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-

# Section D



<b>MODEL:</b> ICX930	
<b>FIELDS:</b> SPOT 3-IN SQ.	
<b>PRE-AMP:</b> 61125B	
<b>OUTPUT:</b> NEGATIVE RAMP	
<b>THICKNESS:</b> 0.4 IN (10.2mm)	
<b>APPLICATION:</b> PICKER 1717 SPOT FILM DEVICE	

<b>MODEL:</b> ICX2822	
<b>FIELDS:</b> A=2, B=1, C=3	
<b>PRE-AMP:</b> 61154S	
<b>OUTPUT:</b> POSITIVE RAMP	
<b>THICKNESS:</b> 0.14 IN (3.5mm)	
<b>APPLICATION:</b> CANON / IDC DIGITAL RADIOGRAPHIC SYSTEM with REMOVABLE REMOTE PRE-AMP	
PICKER COMPATIBLE	

-SERVING THE MEDICAL X-RAY FIELD SINCE 1969-